
How to attach an asset tag to equipment:

- Go to Assets or Tags page in My Site. Locate the asset or tag and select the Attach icon in the Actions column.
- Input the asset and/or asset tag ID in the Tag Attachment window (asset or tag ID field may auto-populate depending on selections in Step 1). Select Next.
- Set the custom properties for the asset.
 Asset details ensure that assets can be easily searched for using a variety of data and identifying information. Click Submit.




Need Assistance? Check out the Cognosos Luminate Healthcare Knowledge Base or contact us at support@cognosos.com.

How to detach an asset tag from equipment:

Go to Assets or Tags in My Site. Locate the asset or tag and choose Detach from the Actions column.



If initiated from Assets, the Asset ID is populated; if from Tags, the Tag ID is populated. Click Confirm to continue.



The tag will have a status of Not Attached. This removes the association between tag and asset; the asset can not be monitored until a new tag is attached.

You may also perform attach/detach in a single step with the Replace Tag option in Asset Details.



Need Assistance? Check out the Cognosos Luminate Healthcare Knowledge Base or contact us at support@cognosos.com.